

# Tech Tips

## Proper Titrating Techniques for Measuring Chemical Concentration



Before working with any chemicals, read and observe precautionary labelling and the MSDS.

Titrations are chemical tests used to determine the concentration of a detergent or sanitizer solution.

Common titrations normally performed in a food plant are:

	<u>Chemical</u>	Reaction	Titrant Indicator	Color Change	
<u>-</u>	Alkali	Acid/Base	Acid (H <sub>2</sub> SO <sub>4</sub> ) Phenolphthalein	Red → Clear	5
-	Acid	Acid/Base	Alkali (NaOH) Phenolphthalein	Clear→ Red	_
	NaOCl or Chlorine	Oxidation - Reduction	Thiosulfate (Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> ) Starch/Iodine	Blue → Clear	
_	Quat (Cationic)	Anionic/ Cationic	Dodecyl SO <sub>4</sub> Chrome Azurol (Anionic)	<b>Purple</b> → <b>Orange</b>	






### **Field Titration Procedure**

Step 1: Collect sample.

Obtain a representative sample from the wash tank, sampling petcock, or line. For obtaining a foam sample, foam the product into a bucket and allow the foam to collapse.

Step 2: Measure Sample. Fill sample vial to calibrated mark. Foam or bubbles do not count. The meniscus is the curved portion of the liquid that should touch the calibrated mark.

Step 3: Titrate. Follow directions in the kit. Hold the bottle straight in a vertical line with the vial and begin titration. Do not hold the bottle at an angle.





### Step 4. Swirl.

Count the drops of the titrant as it is added to the vial. The endpoint is the first permanent color change obtained. A phenolphthalein endpoint should remain permanent for 30 seconds. If needed, repeat titration to confirm results.

Step 5. Calculate results.

The usual formula is:

# Of Drops X Factor of Product = Concentration of Product

Step 6. Record.

Determine if corrective action is required.

Keep in mind, the accuracy of a field test titration is +/- 10%.

#### **Johnson Diversey**

#3000388

3630 East Kemper Road Sharonville, Ohio 45249 1-800-626-5015 2401 Bristol Circle Oakville Ontario, L6H 6PI 1-800-668-7171

